

SaniBact

Medical-surgical disinfectant device, registration N20145

SANIBACT is a disinfectant, which has been designed in Errecom labs for the evaporators, air channels and every kind of surface that need a deep disinfecting level cleaning.

Sanibact works even as a strong cleanser and guarantees the highest level of hygiene during the cleaning process.

Sanibact is a HACCP product and it is ideal for the following applications:

hospitals, nursing home, medical centers, clinic, schools, sport courts, wellness center, butcher shops, sweet companies, bakeries, restaurants, ice-cream shops, hotels, companies that deals with meat, fish and vegetables.

This product does not release odor and/or remains.

Its disinfecting strength acts in just five minutes, allowing a deep and rapid action.

Sanibact is an extremely concentrated product: a dilution of 1:25 is sufficient.

It is efficient against the following virus:

influence, hepatitis B/C, HIV, Rota virus, Noro virus and Ebola Virus.

Concerning bacteria, Sanibact is in compliance with EN1276 EN 13697 tests and for the moulds EN13697 test..

Instruction: Dilute SaniBact 1:25 in water. Spray on the surface with a trigger or with nozzle pump. Rinse after 15 minutes.

Professional usage product: wear the appropriate protective uniform.

Read the instructions on the label carefully.





EN ES

(DE)

(RU)

(cz)

BG

(RO) (GR) (BR)